

Wellbeing Fold Down Toilet Grab Rail 850mm

ECFOLDDOWN

Fold-up rail for use in toilet applications for disabled, aged care facilities and home assist use.

32mm 304 grade stainless steel tubing.



Designed for maximum strength and user confidence.

FEATURES

- Wall mounted
- Folds up when not in use
- 304 grade Satin stainless steel tubing
- 32mm diameter rail 850mm x 210mm

Product Codes

ECFOLDDOWN Fold Down Toilet Grab Rail -
Stainless Steel Finish



Fixing screws are not supplied with this product. Ensure grab rail is installed on a solid wall surface with firm wall backing using appropriate fixing method. Structural engineers or licensed tradespersons should be consulted for appropriate advice on suitability of the fasteners and wall structure before affixing or installing any Enware products. Enware does not provide any warranty in relation to bolts, screws or other fasteners used with Enware products.

NOTE: Enware Australia advises:

1. Due to ongoing Research and Development, specifications may change without notice.
2. Component specifications may change on some export models.
3. Refer to Warranty Statement for warranty details - www.enware.com.au/warranty

Version 1.1__05Sep22

Call 1300 369 273
www.enware.com.au

Enware Australia Pty Limited
9 Endeavour Rd Caringbah NSW 2229 Australia
Ph: 02 8536 4000 info@enware.com.au

