











Ret til ændringer forbeholdes. Recht auf Änderungen sind vorbehalten. Subject to alteration without further notice.
Nous nous réservons le droit d'apporter toutes modifications pour l'amélioration de nos produits sans préavis.
Wijzigingen voorbehouden. Salvo variazioni. Rätten till ändringar förbehålles. Vi förbehåller oss rätt till ändringar.
Kaikki oikeudet muutoksiin pidätetään. Rettur til breytinga áskilinn. Reservado o direito de fazer modificações sem prévio
Επιτηρασόμεθα για εκπαιτωικά γάθη, μικές τεχνικές αλλαγές και παρεκκλίσεις μοντέλων



7a



DK:

Det er vigtigt at bufferne ligger tæt op af toiletskålens inderside.

SE:

Det är viktigt att stabilisera buffertarna så att dessa passar tätt mot insidan av toalett skålen.

EN:

It is important that the stabilizing buffers fit tight against the inner rim of the toilet bowl.


DE:

Es ist wichtig, dass die Puffer eng an der Innenseite der Sanitärkeramik anliegen.

FR:

Il est important qu'il n'y ait aucun jeu entre les plots stabilisateurs et le bord intérieur de la cuvette.

7b



DK:

Det er vigtigt at bufferne ligger tæt op af toiletskålens yderside.

SE:

Det är viktigt att stabilisera buffertarna så att dessa passar tätt mot utsidan av toalett skålen.

EN:

It is important that the stabilizing buffers fit tight against the outside rim of the toilet bowl.


DE:

Es ist wichtig, dass die Puffer eng an der Außenseite der Sanitärkeramik anliegen.


FR:

Il est important qu'il n'y ait aucun jeu entre les plots stabilisateurs et le bord extérieur de la cuvette.


8



9




10




CERA
(DK)



Product	Toilet seat Dania without cover, 50 mm raised		
Item number	R23000 VVS:615684000 NRF:6190526 EAN:5708590070476 HMI: HMI 1811 NBS:N13:300		
Colour	White		
Product description	Toilet seat, Dania, with stabilising buffers, w/o cover, raised 50 mm With raised cross bar hinge		
Adjustment height			
Maximum load	This product has been tested statically up to 225 kg. As the European standard for technical aids for disabled, BS/EN 12182:2012, requires a static test with safety factor 1.5, the nominal maximum load for this product is therefore 150 kg.		
Materials	Seat: Thermoset plastic (urea) Strap: Stainless steel V2A Raiser: Nylon (PA) Special buffer: ABS-plastic		
Ergonomy	An ergonomically designed seat that is well suited to most people. Sideways transfer is easy because the seat is rather flat. The cover is comfortable to sit on, and when set upright it gives good support to the back. Raised 5 cm to make it less difficult to sit down and stand up.		

Product	Toilet seat Dania without cover, 50 mm raised
---------	---



Product	Toilet seat Dania with cover, 50 mm raised	
Item number	R33000 VVS:615685000 NRF:6190584 RSK:7972415 EAN:5708590070520 HMI: HMI 1810 NBS:N13:300	
Colour	White	
Product description	Toilet seat, Dania, with stabilising buffers, w/cover, raised 50 mm With raised cross bar hinge	
Adjustment height		
Maximum load	This product has been tested statically up to 225 kg. As the European standard for technical aids for disabled, BS/EN 12182:2012, requires a static test with safety factor 1.5, the nominal maximum load for this product is therefore 150 kg.	
Materials	Seat/cover: Thermoset plastic (urea) Strap: Stainless steel V2A Raiser: Nylon (PA) Special buffer: ABS-plastic	
Ergonomy	An ergonomically designed seat that is well suited to most people. Sideways transfer is easy because the seat is rather flat. The cover is comfortable to sit on, and when set upright it gives good support to the back. Raised 5 cm to make it less difficult to sit down and stand up.	
Cleaning	Use a mild soap solution to clean the seat. Please take care that the seat and hinges are not left damp. Dry off any residual water or cleaning agent with a soft cloth. Avoid contact of the seat and hinges with abrasive, corrosive or chlorine-based cleaners, as these may cause damage or lead to flash rust. Therefore, when cleaning the bowl, make sure that the seat and cover are kept in an upright position until all cleaning agent has been flushed away.	

Product	Toilet seat Dania with cover, 50 mm raised
---------	--

